A Raman Analysis of Rock Sample for Bull Trout Jade, LLC

Lanny R. Ream Coeur d'Alene Idaho

October 3, 2014

Method and Equipment

Sample was analyzed using a Raman spectrometer with a 532 μ m (green light) laser. The spectra were acquired at 1 second acquisition times, generally taken in a series of 10 acquisitions averaged together. The sample was a small block of blue to blue-gray color with a white area in the center.

Resulting spectra were compared to standards that are available from several sources at universities in the USA and Europe. At this time, about 65% of the known mineral species are represented by at least one Raman standard, but all mineral species normally found in the nephrite deposits are represented in these standards.

Source of Samples

The samples were provided by Rodney Cook of Bull Trout, LLC.

Lanny R. Ream PO Box 2043 Coeur d'Alene, ID 83816

Samples and results

BT8-1

Block of blue jade with variegated white zone in center.

Analyses

BT8-1a bright white zone = calcite

1b creamy white = calcite

1c translucent gray veinlet = calcite

1d white layer = calcite plus titanite

1e white to grayish glob = calcite plus titanite


1f translucent gray zone = calcite

1g medium blue gray = tremolite


1h light blue gray = tremolite

1i light gray = titanite


1j light gray in calcite mass = titanite plus calcite




Analysis map of sample BT8-1.



BT8-1a spectrum, black is sample, blue is a calcite standard



BT8-1d spectrum, black is sample, blue is a calcite standard green is a titanite standard



BT8-1g spectrum, black is sample, blue is a tremolite standard